reprocessing of Polydentia products

Disinfection,
Cleaning &
Sterilisation
instructions

polydentia
swiss manufacture 

Index

EN - English	2
DE - Deutsch	6
FR - Français	11
IT - Italiano	16
ES - Español	21

Reprocessing of Polydentia Products

Disinfection, Cleaning & Sterilisation instructions

Basic Principles:

- All Polydentia reusable Products are supplied non-sterile and must be appropriately disinfected, cleaned and sterilized before the first use.
- All Polydentia Products must to be disinfected, cleaned and sterilized prior to each use.
- Effective cleaning is an indispensable requirement for proper Polydentia Products sterilization.
- Consider the legal provisions valid for your country as well as to the hygienic instructions of the Dental Practice.
- The user is responsible for the sterility of Polydentia Products. Therefore, please ensure that only validated procedures are used for disinfection, cleaning and sterilization. The sterilization equipment must also be maintained and checked regularly, as well as the validated parameters applied to each disinfection, cleaning and sterilization cycle.
- Materials may alter over time. Sterilization or exposure to chemicals can accelerate this deterioration. Always check your Polydentia Product for wear and damage before use.
- Avoid contact of dissimilar metals at any point during the disinfection, cleaning and sterilization process.

Protection of Staff Members:

All used and contaminated Polydentia Products must be handled with appropriate Personal Protective Equipment.

1. Pre-treatment

Before processing Polydentia Products remove coarse impurities on the products immediately after application.

Do not use a metal brush. A metal brush will cause scratches or damages on the products surface.

- ▶ In case of using a thermal disinfector the steps in this guidance (see below) “2. Disinfection” and “3. Rinsing” are unnecessary, provided that the dirty Polydentia Products have not been stored dry longer than 2 hours.

2. Pre-treatment

- Do not leave dirty Polydentia Products dry – put them into a disinfecting liquid as soon as possible.
- We suggest to use an enzymatic pre-cleaner disinfectant such as, but not limited to:
 - **Dürr Dental ID215**
 - Or a pre-cleaning disinfectants based on Quaternary Ammonium compounds.
- Follow carefully the instructions for use given by the manufacturer of the selected disinfectant. Please note that too strong concentration of the disinfectant or too long exposure time may damage the products.
- Change the disinfectant regularly according to the instructions.

3. Rinsing

- Immediately after the disinfection, Polydentia Products must be rinsed carefully under running lukewarm water.

4. Cleaning

Always clean Polydentia Products, disinfection and rinsing are not sufficient.

Cleaning alternatives:

Manual (A), Ultrasonic bath (B) Thermal Disinfector (C)

If possible, an automatic procedure in a thermal disinfector should be used for cleaning of Polydentia Products. A manual procedure, even in case of application of an ultrasonic bath, should only be used if an automatic procedure is not available or if such a method is not compatible with specific materials; in this case, the significantly lower efficiency of a manual procedure must be considered.

A. Manual Cleaning

- **Completely disassemble Polydentia Products if applicable**
- We suggest to use a cleaning solution such as, but not limited to:
 - **Dürr Dental ID212**
 - Or cleaning solutions based on Quaternary Ammonium compounds. Check the pH-value of the cleaning solutions; it must be as neutral as possible.
- Follow carefully the instructions for use given by the manufacturer of the selected cleaning solutions. Please note that too strong concentration or too long exposure time may damage Polydentia Products.
- Change the cleaning solutions regularly according to the instructions.
- Always rinse carefully Polydentia Products under running water after cleaning.

B. Ultrasonic bath

- Please assure that Polydentia Products have the relevant symbol on the packaging.
- Completely disassemble Polydentia Products if applicable.
- Follow carefully the instructions for use of the ultrasonic bath device.
- Service and clean the device frequently, according to the instructions for use.
- Ensure that maximum load is not exceeded.
- We recommended to use a cleaning solutions such as, but not limited to:
 - Dürr Dental ID212
 - Or cleaning solutions based on Quaternary Ammonium compounds. Check the pH-value of the cleaning solutions; it must be as neutral as possible.
- Follow carefully the instructions for use given by the manufacturer of the selected cleaning solutions. Please note that too strong concentration or too long exposure time may damage Polydentia Products.
- Change the cleaning solutions regularly according to the instructions.
- Remove the devices from the ultrasonic bath immediately after cleaning.
- Always rinse carefully Polydentia Products under running water after cleaning.

C. Thermal Disinfector

- Please assure that Polydentia Products have the relevant symbol on the packaging.
- Completely disassemble Polydentia Products if applicable.
- Follow carefully the instructions for use given by the manufacturer of the thermal disinfector and of the selected program; and prefer cleaning agents with corrosion protector.
- Inspect the device regularly according to the instructions for use. Service and clean the device frequently, according to the instructions for use.
- Ensure that maximum load is not exceeded.
- Ensure that the selected washing programme and cleaning agents are effective enough.
- Remove Polydentia Products from the thermal disinfector immediately after cleaning.

5. Drying

Before sterilizing in the autoclave, Polydentia Products must be thoroughly dried. The normal tap water may contain different chemical ingredients which may cause chemical reactions when autoclaving.

6. Checking

After cleaning, check that Polydentia Products are clean by rubbing the surface with an eraser or finger.

WARNING! It is extremely important to check that the products are truly clean and dry before the sterilizing procedure.

7. Sterilization

Sterilize Polydentia Products **only** in a steam autoclave using distilled water and at the temperature specified on the relevant symbol. **Do not use chemical or cold sterilization.**

- **Please assure that Polydentia Products have the relevant symbol on the packaging.**
- Follow carefully the instructions for use of the autoclave.
- Inspect the device regularly according to the instructions for use. Service and clean the device frequently, according to the instructions for use .
- Ensure that maximum load is not exceeded.

- Place all the medical device components in an autoclave tray, or in a sterilization disposable pouch (recommendet).
- Plastic parts must not touch the autoclave walls, since the temperature there might be higher.
- We recommended to use always the following cycle (Prion Cycle):
Temperature of sterilisation: 134°C
Pressure of sterilisation: 2 bar
Time of sterilisation: 20.5 min
Drying time: 20 min

- Remove Polydentia Products from the autoclave immediately after sterilisation.

8. Storage

- Store sterilized product in a dry and dust free location.
- Keep product sealed in the sterilization pouch until ready to use.

DISCLAIMER:

This high quality product is manufactured by Polydentia SA, in Mezzovico, Switzerland. The content of this packaging has been inspected, tested and approved by our quality officers. Materials used during our manufacturing process are designed to provide durability and performance. Polydentia constantly looks for the highest quality plastics for its product; nevertheless please take into account that plastics may alter over time. Sterilization or exposure to chemicals can accelerate this deterioration. For the longest lifetime of your premium quality product always clean and sterilize it following these instructions.

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