Instructions for use



F-Splint-Aid F-Splint-Aid Slim

splinting



Index

EN	English	6
DE	Deutsch	8
FR	Français	10
IT	Italiano	12
ES	Español	14
PT	Português	16
NL	Nederlands	18
NO	Norsk	20
FI	Suomi	22
SE	Svenska	24
DK	Dansk	26
PL	Polski	28
CZ	Čeština	30
TK	Türkçe	32
GR	ελληνικά	34
JP	日本語	36
CN	中国	38
AR	عربي	40



I






II

III

IV

A. DESCRIPTION

- **F-Splint-Aid** Fiberglass strip (6 layers, 4mm wide), impregnated with Fiber-Bond bonding agent.
- **F-Splint-Aid Slim** Fiberglass strip (3 layers, 2mm wide), impregnated with Fiber-Bond bonding agent.

• Bonding component Fiber-Bond:

- without Primer;
- not radio-opaque;
- suitable for the "Etching the enamel - Etch And Rinse" bonding technique;

B. COMPOSITION

Strip Fiberglass

Bonding agent (Fiber-Bond)

2,2-Bis-4 (2hydroxide-3-methacryloxyde-propyloxyde) - phenol - propane

2,2-Bis-4 (2-methacryloxyde-ethoxyde) - phenol - propane

3,6-dioxy octamethylene dimethacrylate

C. INDICATIONS

- Stabilization of teeth with increased mobility due to periodontitis.
- Stabilization of teeth with increased mobility or avulsed following a trauma.
- Adhesive front bridges created with teeth in composite build-up material, with prefabricated teeth and original teeth with resected roots.

D. CONTRAINDICATIONS

Do not use F-Splint-Aid/F-Splint-Aid Slim on patients known to be allergic to one of the components.

E. SIDE EFFECTS

For the intended use in which F-Splint-Aid/F-Splint-Aid Slim do not come in contact with the dentin, no side effects have yet been observed.

F. INTERACTIONS WITH OTHER DEVICES

- F-Splint-Aid/F-Splint-Aid Slim must not be used with products containing eugenol since these can inhibit the curing process.
- Do not use any resin to modify the viscosity of the bonding component of F-Splint-Aid/F-Splint-Aid Slim.
- For the correct positioning of F-Splint-Aid/F-Splint-Aid Slim on the teeth, use Clip&Splint application clips.

G. WARNINGS

Since F-Splint-Aid/F-Splint-Aid Slim are sensitive to blue light, curing is activated also by daylight or artificial light (lighting of the operating field). The curing process can be slowed down considerably by using filters that inhibit the passage of blue light. It is therefore advisable to expose F-Splint-Aid/F-Splint-Aid Slim to light only for the time strictly necessary (pick up the fiber strip only at the time of application and close the bottle immediately after taking the strip) and apply it rapidly.

Incorrect application procedure increases the risk of device detachment.

H. APPLICATION

See pictograms

I. Preparing the teeth for application of the Fiber-Bond bonding agent

1. After scaling, clean all the affected teeth thoroughly with a prophylactic non-fluoride paste and apply a dam.
2. Etch the palatal/lingual and interdental surface with 37% orthophosphoric acid according to the instructions for use provided by the manufacturer of the etching product used. Existing fillings and/or restorations can be included in the preparation if etched appropriately.

II. Applying Fiber-Bond

3. If necessary, highly mobile teeth can be stabilized temporarily by applying a few drops of Fiber-Bond or fluid composite (by following the instructions for use provided by the manufacturer of the composite) in the interdental spaces on the buccal side.

- Cover the palatal/lingual and interdental surfaces to be treated with a layer of Fiber-Bond. A few drops of bonding agent can be taken directly from the F-Splint-Aid/F-Splint-Aid Slim bottle by pressing it slightly.

III. Preparing the fiberglass strip

- Using a pair of tweezers take a strip of suitable length from the bottle of F-Splint-Aid/F-Splint-Aid Slim.
- Position the bottle with the fiber facing downwards and press slightly on the sides of the bottle to further saturate the fiberglass strip with Fiber-Bond.
- Cut the impregnated strip to the required length, taking care to leave at least 2 mm of strip coming out of the nozzle to make subsequent use easier. Close the bottle of F-Splint-Aid/F-Splint-Aid Slim perfectly after use to preserve the product.

IV. Positioning the fiberglass strip on the teeth, light curing and finishing

- Position the strip on the palatal/lingual surface and hold it in place by means of the Clip&Splint application clips applied in the interdental spaces.
- Hold the ends of the strip in place by means of tweezers and cure each interdental space for the time indicated below.
Hardening of the fiberglass strip (when using the bonding agent Fiber-Bond) depends on the performance of the curing lamp. As a rule, for a 1 cm long strip a curing time of 20 seconds is sufficient. However, to ensure uniformity, approx. 30-40 seconds is recommended. This estimate is based on assumption of a light intensity of approx. 800 mW/cm², which corresponds to the usual intensity of lamps having wavelength of 400-500 nm (blue light). When curing with LED light, since many LED units emit light with wavelength of only 430 nm, a slightly higher power is recommended.
- Any irregularities or roughness formed must be eliminated by applying a thin layer of fluid composite, according to the instructions for use provided by the manufacturer of the composite used. Refine and polish. Remove the dam and check the occlusion carefully.

I. FIRST AID MEASURES

- In case of contact with the eyes rinse immediately and thoroughly with running water for 10-15 minutes and consult an ophthalmologist if necessary.
- In case of contact with the skin, rinse immediately and thoroughly. See a doctor if necessary.

J. INFORMATION FOR THE PATIENT

It is very important to inform the patient of the correct prophylaxis to be adopted after the application of F-Splint-Aid/F-Splint-Aid Slim.

K. SHELF LIFE

F-Splint-Aid/F-Splint-Aid Slim must not be used beyond the expiration date indicated on the package and on the bottle.

L. PRESERVATION AND STORAGE

Storage in the refrigerator at +5°C (41°F) to +8°C (47°F) is recommended.

Avoid temperatures higher than 25°C (77°F).

M. PACK SIZES AVAILABLE

See the "Reorder now" Table

N. DISPOSAL

Disposal of F-Splint-Aid/F-Splint-Aid Slim must be done in accordance with the applicable laws and regulations.

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reorder now

REF 5972	F-Splint-Aid: Fiber-Splint ribbon (4 mm wide, 120 mm long) impregnated with Fiber-Bond + 5 Clip&Splint application clips
REF 5977	F-Splint-Aid Slim: Fiber-Splint ribbon (2 mm wide, 120 mm long) impregnated with Fiber-Bond + 5 Clip&Splint application clips
REF 5980	F-Splint-Aid & Slim Introkit: REF 5972 + REF 5977
REF 5975CS	Clip&Splint, autoclavable application clips, 30 pcs
REF 6525	Fiber-Bond: 2 bottle x 5ml

Symbols

REF	Catalogue number
LOT	Batch code
	Use by date
	Consult instructions for use
	Do not reuse
 -25 °C (77 °F)	Upper limit of temperature
Rx Only	For use only by trained medical personnel
	Manufacturer
CE 1639	CE - Notified body identification number

Importer

Unique device identifier

Authorized representative in the
European Community/ European Union

Rx Only CE 1639

Only for professional use
REV. 06.01.2022

 Polydentia SA
via Cantonale 47
6805 Mezzovico-Vira
Switzerland

Made in Switzerland

phone: +41 91 946 29 48
fax: +41 91 946 32 03

info@polydentia.ch
www.polydentia.ch