Hemo- Wedges

Instructions for use

Index

EN - English	2
DE - Deutsch	3
FR - Français	4
IT - Italiano	5
ES - Español	6

HEMO-WEDGES

Interdental wedges of maple wood with hemostatic properties

COMPOSITION

Maple wood, impregnated with aluminium potassium sulphate in the following concentration:

**Mass Percentage Rate [w%] =
7% of the total weight of each
wedge**

PROPERTIES

The aluminium potassium sulphate, which is contained in the wedge, has an astringent and hemostatic effect on the surface of the gingiva. The effect covers the whole period in which the interdental wedge and the gingiva remain in contact, either in a direct or in a fluid-induced way.

INDICATIONS

All cases requiring a conservative treatment with composite, presenting serous and bloody exudations.

APPLICATION MODALITIES

The application modalities of the device are identical to those of the usual interdental wood wedges.

CONTRAINDICATIONS

Eventual allergies to aluminium potassium sulphate.

WARNINGS AND SAFETY PRECAUTIONS

- Thoroughly rinse the site of the contact with the interdental wedge after its removal. Eventual residuals of aluminium potassium sulphate may cause irritation.
- Avoid hand-eye contact after using the device, because this may have irritant effect.

STORAGE

- Store at room temperature (15-25 °C), in a dry place, protected from light and from heat.
- Do not use the device beyond the expiration date indicated on the packaging

PRESENTATION

REF list on the product label.

PRESCRIPTION CONDITIONS

For professional dental use only.

REVISION

See revision index on the present instructions.

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