Instructions for use

Fiber-Splint  
Fiber-Splint ML Multi-Layer  
Fiber-Splint Ortho Evolution  
Fiber-Bond

splinting

polydentia  
swiss manufacture 

# index



EN	English	6
DE	Deutsch	12
FR	Français	18
IT	Italiano	24
ES	Español	30
PT	Português	36
NL	Nederlands	42
NO	Norsk	48
FI	Suomi	54
SE	Svenska	60
DK	Dansk	66
PL	Polski	72
CZ	Čeština	78
TK	Türkçe	84
GR	ελληνικά	90
JP	日本語	96
CN	中国	102
AR	عربي	107


I


II


III

IV

V

### IMPORTANT INFORMATION

When applying Fiber-Splint/Fiber-Splint ML Multi-Layer/Fiber-Splint Ortho Evolution according to paragraph H of these instructions, the use of bonding agent Fiber-Bond by Polydentia is recommended.

The Fiber-Bond bonding agent must be used only and exclusively in combination with the products, Fiber-Splint/Fiber-Splint ML Multi-Layer/Fiber-Splint Ortho Evolution for the intended use described in paragraph C of these instructions, where Fiber-Bond is only applied on the enamel. Fiber-Bond must not be used for any other purpose.

#### A. DESCRIPTION

##### Fiber-Splint

Fiberglass strip (single layer, 4mm wide)

##### Fiber-Splint ML Multi-Layer

Fiberglass strip (6 layers, 4mm wide)

##### Fiber-Splint Ortho Evolution

Fiberglass strip (3 layers, 2mm wide)

##### Fiber-Bond

Bonding agent:

- without Primer;
- not radio-opaque;
- suitable for the "Etching the enamel - Etch And Rinse" bonding technique;

#### B. COMPOSITION

Fiber-Splint	Fiberglass
Fiber-Splint ML Multi-Layer	Fiberglass
Fiber-Splint Ortho Evolution	Fiberglass

##### Fiber-Bond (bonding agent)

2,2-Bis-4 (2hydroxide-3-methacryloxyde-propyloxyde)  
- phenol - propane

2,2-Bis-4 (2-methacryloxy-ethoxy) - phenol -  
propane  
3,6-dioxy octamethylene dimethacrylate

### C. INDICATION

- Stabilization of teeth with increased mobility due to periodontitis.
- Stabilization of teeth with increased mobility or avulsed following a trauma.
- Adhesive front bridges created with teeth restored in composite, with prefabricated teeth and original teeth with resected roots.

### D. CONTRAINDICATIONS

Do not use Fiber-Splint/Fiber-Splint ML Multi-Layer/  
Fiber-Splint Ortho Evolution/Fiber-Bond on patients  
known to be allergic to one of the components.

### E. SIDE EFFECT

For the intended use in which Fiber-Splint/Fiber-Splint ML Multi-Layer/Fiber-Splint Ortho Evolution in combination with Fiber-Bond do not come in contact with the dentin, no side effects have yet been observed.

### F. INTERACTIONS WITH OTHER DEVICES

- The Fiber-Bond bonding agent must not be used together with products containing eugenol since these can inhibit the curing process.
- Do not use any resin to modify the viscosity of the Fiber-Bond bonding agent.
- To position Fiber-Splint/Fiber-Splint ML Multi-Layer/Fiber-Splint Ortho Evolution correctly on teeth, use Clip&Splint application clips.

### G. WARNINGS

Since the Fiber-Bond bonding agent is sensitive to blue light, curing is activated also by daylight or artificial

light (lighting of the operating field). The curing process can be slowed down considerably by using filters that inhibit the passage of blue light. It is therefore advisable to expose the Fiber-Bond bonding agent to light only for the time strictly necessary for the preparation and application of Fiber-Splint/Fiber-Splint ML Multi-Layer/Fiber-Splint Ortho Evolution. The Fiber-Bond bottle must be closed every time after use.

Incorrect application procedure increases the risk of device detachment.

## H. APPLICATION

See pictograms

### I. Preparing the teeth for application of the Fiber-Bond bonding agent

1. After scaling, clean all the affected teeth thoroughly with a prophylactic non-fluoride paste and apply a dam.

2. Etch the palatal/lingual and interdental surface with 37% orthophosphoric acid according to the instructions for use provided by the manufacturer of the etching product used. Existing fillings and/or restorations can be included in the preparation if etched appropriately.

### II. Applying the bonding agent

**After etching the enamel it is not necessary to use a Primer before applying Fiber-Bond.**

3. If necessary, highly mobile teeth can be stabilized temporarily by applying a few drops of Fiber-Bond or another bonding agent or fluid composite (by following the instructions for use provided by the manufacturer of the composite) in the interdental spaces on the buccal side.
4. Cover the palatal/lingual and interdental surfaces to be treated with a layer of Fiber-Bond or another bonding agent.

### **III. Preparing the fiberglass strip**

#### **5a. Fiber-Splint**

Cut 6 strips to the required length.

#### **5b. Fiber-Splint ML Multi-Layer / Fiber-Splint Ortho Evolution**

Cut the strip to the required length.

### **IV. Impregnating the fiberglass strip with bonding agent**

#### **6a. Fiber-Splint**

Impregnate the first strip of Fiber-Bond or other bonding agent on a glass plate. The fabric must be completely soaked.

#### **6b. Fiber-Splint ML Multi-Layer / Fiber-Splint Ortho Evolution**

Impregnate the strip of Fiber-Bond or other bonding agent on a glass plate. The fabric must be completely soaked.

### **V. Positioning the fiberglass strips on the teeth, light curing and finishing**

#### **7a. Fiber-Splint**

Position the first strip on the palatal/lingual surface and hold it in place by means of the Clip&Splint application clips applied in the interdental spaces.

#### **7b. Fiber-Splint ML Multi-Layer / Fiber-Splint Ortho Evolution**

Position the strip on the palatal/lingual surface and hold it in place by means of the Clip&Splint application clips applied in the interdental spaces.

#### **8. Hold the ends of the strip in place by means of tweezers and cure each interdental space for the time indicated below.**

Hardening of the fiber strip applied with the Fiber-Bond bonding agent depends on the characteristics of the curing lamp. As a rule, for a 1 cm long strip a curing time of 20 seconds is

sufficient. However, to ensure uniformity, approx. 30-40 seconds is recommended. This estimate is based on a light intensity of about 800 mW/cm<sup>2</sup>, which corresponds to the usual intensity of lamps with wavelength 400-500 nm (blue light). When curing with LED light, since many LED units emit light with wavelength of only 430 nm, a slightly higher power is recommended.

### **Fiber-Splint**

Repeat steps 6 to 9 for the remaining 5 strips.

9-11. Any irregularities or roughness must be eliminated formed by applying a thin layer of fluid composite, according to the instructions for use provided by the manufacturer of the composite used. Refine and polish. Remove the dam and check the occlusion carefully.

### **I. FIRST AID MEASURES**

- If the Fiber-Bond bonding agent comes in contact with the eyes, rinse immediately and thoroughly with plenty of running water for 10-15 minutes and consult an ophthalmologist if necessary.
- If the Fiber-Bond bonding agent comes in contact with the skin, rinse thoroughly immediately. See a doctor if necessary.

### **J. INFORMATION FOR THE PATIENT**

It is very important to inform the patient of the correct prophylaxis to be adopted after the application of Fiber-Splint/Fiber-Splint ML Multi-Layer/Fiber-Splint Ortho Evolution.

### **K. SHELF LIFE**

Fiber-Splint/Fiber-Splint ML Multi-Layer/Fiber-Splint Ortho Evolution/Fiber-Bond must not be used beyond

the expiration date indicated on the package and on the bottle.

#### L. PRESERVATION AND STORAGE

Only for the Fiber-Bond bonding agent storage in a refrigerator from +5°C (41°F) to +8°C (47°F) is recommended.

Avoid temperatures higher than 25°C (77°F).

#### M. PACK SIZES AVAILABLE

See the “Reorder now” Table.

#### N. DISPOSAL

Disposal of Fiber-Splint/Fiber-Splint ML Multi-Layer/Fiber-Splint Ortho Evolution/Fiber-Bond must be done in accordance with the applicable laws and regulations.

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## reorder now

REF 5970	Fiber-Splint (4 mm wide, 200 cm long) + 5 Clip&Splint application clips
REF 5971	Fiber-Splint ML Multi-Layer (4 mm wide, 50 cm long) + 5 Clip&Splint application clips
REF 5979	Fiber-Splint Ortho Evolution (2 mm wide, 60 cm long) + 5 Clip&Splint application clips
REF 5975CS	Clip&Splint, autoclavable application clips, 30 pcs
REF 6525	Fiber-Bond: 2 bottle x 5ml
REF 5980	F-Splint-Aid & Slim Introkit: REF 5972 + REF 5977
REF 5972	F-Splint-Aid (4 mm wide, 120 cm long) + 5 Clip&Splint application clips
REF 5977	F-Splint-Aid Slim (2 mm wide, 120 cm long) + 5 Clip&Splint application clips

## Symbols

<b>REF</b>	Catalogue number
<b>LOT</b>	Batch code
	Use by date
	Consult instructions for use
	Do not reuse
 <sup>25 °C (77 °F)</sup>	Upper limit of temperature
<b>R<sub>x</sub> Only</b>	For use only by trained medical personnel
	Manufacturer
<b>MD</b>	Medical Device
<b>CE</b> 1639	CE - Notified body identification number

	Importer
<b>UDI</b>	Unique device identifier
<b>EC</b> <b>REP</b>	Authorized representative in the European Community / European Union

Only for Fiber-Bond

R<sub>X</sub> Only

MD

CE 1639

Only for professional use

REV. 07.04.2022

**polydentia**  
swiss manufacture 

Polydentia SA  
via Cantonale 47  
6805 Mezzovico-Vira  
Switzerland

Made in Switzerland

phone: +41 91 946 29 48  
fax: +41 91 946 32 03

info@polydentia.ch  
www.polydentia.ch